

# Establishment Registration & Device Listing



FDA Home | Medical Devices | Databases

New Search		Back To Search Results
<b>Proprietary Name:</b>	AIR QUEEN; AIR QUEEN BREEZE; Nano Mask; PURE MASK; TECHNO WEB; TECHNWEB	
<b>Classification Name:</b>	MASK, SURGICAL	
<b>Product Code:</b>	<a href="#">FXX</a>	
<b>Device Class:</b>	2	
<b>Regulation Number:</b>	<a href="#">878.4040</a>	
<b>Medical Specialty:</b>	General & Plastic Surgery	
<b>Registered Establishment Name:</b>	<a href="#">NANO FILTER, INC</a>	
<b>Registered Establishment Number:</b>	3016766213	
<b>Premarket Submission Number:</b>	<a href="#">K172500</a>	
<b>Owner/Operator:</b>	<a href="#">TOPTEC Co., LTD.</a>	
<b>Owner/Operator Number:</b>	10072713	
<b>Establishment Operations:</b>	Contract Manufacturer	

Page Last Updated: 06/08/2020

Note: If you need help accessing information in different file formats, see [Instructions for Downloading Viewers and Players](#).

Language Assistance Available: [Español](#) | [繁體中文](#) | [Tiếng Việt](#) | [한국어](#) | [Tagalog](#) | [Русский](#) | [العربية](#) | [Kreyòl Ayisyen](#) | [Français](#) | [Polski](#) | [Português](#) | [Italiano](#) | [Deutsch](#) | [日本語](#) | [فارسی](#) | [English](#)

## U.S. Food and Drug Administration

10903 New Hampshire Avenue  
Silver Spring, MD 20993  
Ph. 1-888-INFO-FDA (1-888-463-6332)  
[Contact FDA](#)





# 510(k) Premarket Notification

◀ FDA Home ▶ Medical Devices ▶ Databases



510(k) | DeNovo | [Registration & Listing](#) | [Adverse Events](#) | [Recalls](#) | [PMA](#) | [HDE](#) | [Classification](#) | [Standards](#)  
[CFR Title 21](#) | [Radiation-Emitting Products](#) | [X-Ray Assembler](#) | [Medsun Reports](#) | [CLIA](#) | [TPLC](#)

[New Search](#)

[Back To Search Results](#)

<b>Device Classification Name</b>	<a href="#">Mask, Surgical</a>
<b>510(K) Number</b>	K172500
<b>Device Name</b>	Technoweb Surgical Mask
<b>Applicant</b>	YTS GLOBAL INC. 7406 ALBAN STATION CT STE A 108 Springfield, VA 22150
<b>Applicant Contact</b>	Eddie Nguyen
<b>Correspondent</b>	YTS GLOBAL INC. 7406 ALBAN STATION COURT SUITE A108 Springfield, VA 22150
<b>Correspondent Contact</b>	Eddie Nguyen
<b>Regulation Number</b>	<a href="#">878.4040</a>
<b>Classification Product Code</b>	<a href="#">FXX</a>
<b>Date Received</b>	08/18/2017
<b>Decision Date</b>	03/01/2018
<b>Decision</b>	Substantially Equivalent (SESE)
<b>Regulation Medical Specialty</b>	General & Plastic Surgery
<b>510k Review Panel</b>	General Hospital
<b>Summary</b>	<a href="#">Summary</a>
<b>Type</b>	Traditional
<b>Reviewed By Third Party</b>	No
<b>Combination Product</b>	No

Page Last Updated: 06/08/2020

Note: If you need help accessing information in different file formats, see [Instructions for Downloading Viewers and Players](#).

Language Assistance Available: [Español](#) | [繁體中文](#) | [Tiếng Việt](#) | [한국어](#) | [Tagalog](#) | [Русский](#) | [العربية](#) | [Kreyòl Ayisyen](#) | [Français](#) | [Polski](#) | [Português](#) | [Italiano](#) | [Deutsch](#) | [日本語](#) | [فارسی](#) | [English](#)



**FDA U.S. FOOD & DRUG  
ADMINISTRATION**

May 18, 2018

Yts Global Inc.  
Eddie Nguyen  
Senior Logistics Specialist  
7406 Alban Station Court Suite A108  
Springfield, Virginia 22150

Re: K172500  
Trade/Device Name: Technoweb Surgical Mask  
Regulation Number: 21 CFR 878.4040  
Regulation Name: Surgical apparel  
Regulatory Class: Class II  
Product Code: FXX  
Dated: January 23, 2018  
Received: January 29, 2018

Dear Eddie Nguyen:

This letter corrects our substantially equivalent letter of March 01, 2018.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

 Tina Kiang  
The signature is a stylized, cursive script of the name "Tina Kiang" written in black ink. To the right of the signature is a large, bold, black graphic of the letters "FDA" in a serif font, with the "Tina Kiang" text overlaid on it.

Tina Kiang, Ph.D.  
Acting Director  
Division of Anesthesiology,  
General Hospital, Respiratory,  
Infection Control, and Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K172500

Device Name  
Technoweb Surgical Mask

### Indications for Use (Describe)

The Technoweb Surgical Mask is intended to be worn by operating room personnel during surgical procedures to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRAStaff@fda.hhs.gov](mailto:PRAStaff@fda.hhs.gov)



*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) SUMMARY

The assigned 510(k) number is : K172500

<b>510(k) Owner/Applicant</b>	Finetex EnE, inc
<b>US Correspondent</b>	Eddie Nguyen YTS GLOBAL INC 7406 ALBAN STATION COURT SUITE A108 SPRINGFIELD, VA 22150 (703)340-8178
<b>Date Prepared</b>	February 26, 2018
<b>Trade Name</b>	Technoweb Surgical Mask,
<b>Common Name</b>	Surgical Respirator
<b>Classification Name</b>	Surgical Apparel
<b>Review Panel</b>	General Hospital
<b>Product Code</b>	MSH
<b>Device Classification</b>	Class II per 21 CFR 878.4040
<b>Predicate Device</b>	3M VFlex™ Health Care Particulate Respirator and Surgical Mask, Model 1805/1805S cleared under K121069

Table – Substantial Equivalence Comparison.

Proprietary		Finetex EnE, inc	3M Health Care	Substantial Equivalence
Common Name		TechnoWeb Surgical Mask	3M™ VFlex™ HealTh Care Particulate Respirator and Surgical mask,	SAME
Model Number		TM-R, TM-S	1805/1805S	-
Manufacturer		Finetex EnE, inc	3M Health Care	-
510(k) Reference Number		K172500	K121069	
Intended use		Surgical Mask/Respirator	Surgical Mask/Respirator	SAME
Design				SIMILAR
Single Use		Yes	Yes	SAME
Material	Outside cover web	Felt	Polypropylene spunbond	SIMILAR
	Stiffener web	N/A	Polypropylene spunbond	SIMILAR
	Filter web	Felt	Polypropylene	SIMILAR
	Inner web	Felt	Polypropylene	SIMILAR
	Nose-Clip	Aluminum	Aluminum	SAME
	Staple	N/A	Steel	SIMILAR
	Headband	Nylon	Polyisoprene	SIMILAR
	Nose foam	Plastic coating	N/A	SIMILAR
Dimensions		TM-R is 212mm(W), 74mm(H) and model TM-S is 177.8 mm(W), 138 mm(H).		SIMILAR
Product style		Flat fold	Flat fold	SAME
Design Features		Ear loops	Advanced Electrostatic Media	SIMILAR
Fluid Resistance (ASTM F1862)		fluid resistance at 12 mmHg.	YES	SIMILAR

Indication for Use	is intended to be worn by operating room personnel during surgical procedures to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.	is intended to be worn by operating room personnel during surgical procedures to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.	SAME
--------------------	---	---	------

**Device Description** Technoweb Surgical Mask is a type of surgical mask covers the user's nose and mouth. And provides a physical barrier to fluids and particulate materials. Technoweb Surgical Mask is a flat, pleated tie-on or elastic ear-loop mask consisting of three nonwoven layers; inner and outer cover web with middle different filter web sandwiched in between.

**Intended Use** The Technoweb Surgical Mask is intended to be worn by operating room personnel during surgical procedures to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.

**Model Numbers** TM-R and TM-S (Different Sizes)

**Technological Characteristics** Technoweb Surgical Mask is a type of surgical mask covers the user's nose and mouth. And provides a physical barrier to fluids and particulate materials. Technoweb Surgical Mask is a flat, pleated tie-on or elastic ear-loop mask consisting of three nonwoven layers; inner and outer cover web with middle different filter web sandwiched in between.

**Performance Testing** Laboratory testing regarding characteristics was performed on Technoweb surgical mask to verify its safety and performance. A biocompatibility assessment was performed on the patient contact materials of Technoweb surgical mask.

**Conclusions** Finetex EnE, inc concludes that the Technoweb surgical mask is as safe and as effective and is substantially equivalent to predicate device, 3M Health Care Particulate Respirator and Surgical Mask,